



# SaeboGlove | Order Form



Questions? Call Lyncare 061-379090 | Email this order form to: [info@lyncare.ie](mailto:info@lyncare.ie)  
 Lyncare Ltd | Unit 3 Shannonside Business Park | Birdhill | Co. Tipperary | V94 V8W3

Name of Purchaser:		Phone:	
Date of Order:		Purchase Order Number :	
<b>SHIPPING INFORMATION</b>			
Hospital/ Clinic Name:		For Attention of:	
Department:		Address:	
City:	County:	Eircode:	
<b>INVOICE INFORMATION</b>			
Hospital/ Clinic Name:		For Attention of:	
Street:		Dept./Suite:	
City:	County:	Eircode:	

**ORDER BELOW: PLEASE REFER TO THE SIZING INSTRUCTIONS ON PAGE 2.  
 IDENTIFY APPROPRIATE SIZE BY DETERMINING BOTH WRIST SPLINT AND GLOVE LINER MEASUREMENTS IN CM.**

Wrist Splint Sizing in CM

Glove Liner Sizing in CM

		(A) WRIST CIRCUMFERENCE					
		10.2-11.4	11.4-14	14-16.5	16.5-19	19-21.6	>21.6
(B) DORSAL HAND WIDTH	<6.3	S	S	S			
	6.3-7.6	S	S	M	M		
	7.6-8.9	S	M	M	M	M	L
	8.9-10.2			M	M	L	L
	>10.2				L	L	L

		(C) HAND CIRCUMFERENCE					
		11.4-14	14-16.5	16.5-19	19-21.6	21.6-24.1	24.1-26.7
(D) 3 <sup>RD</sup> FINGER LENGTH	15.2-16.5	XS	XS	XS			
	16.5-17.8	XS	XS	XS	S		
	17.8-19		S	S	S	M	
	19-20.3			M	M	L	XL
	>20.3				L	XL	XL

	WRIST SPLINT SIZE	GLOVE LINER SIZE	PRODUCT CODE	QTY	PRICE	TOTAL	
LEFT	Small	Extra Small	G1-GL-SM-XS-L-1		€399.00		
	Small	Small	G1-GL-SM-SM-L-1		€399.00		
	Medium	Extra Small	G1-GL-MD-XS-L-1		€399.00		
	Medium	Small	G1-GL-MD-SM-L-1		€399.00		
	Medium	Medium	G1-GL-MD-MD-L-1		€399.00		
	Medium	Large	G1-GL-MD-LG-L-1		€399.00		
	Large	Large	G1-GL-LG-LG-L-1		€399.00		
	Large	Extra Large	G1-GL-LG-XL-L-1		€399.00		
RIGHT	Small	Extra Small	G1-GL-SM-XS-R-1		€399.00		
	Small	Small	G1-GL-SM-SM-R-1		€399.00		
	Medium	Extra Small	G1-GL-MD-XS-R-1		€399.00		
	Medium	Small	G1-GL-MD-SM-R-1		€399.00		
	Medium	Medium	G1-GL-MD-MD-R-1		€399.00		
	Medium	Large	G1-GL-MD-LG-R-1		€399.00		
	Large	Large	G1-GL-LG-LG-R-1		€399.00		
	Large	Extra Large	G1-GL-LG-XL-R-1		€399.00		
<p style="color: red;"><b>Price includes VAT @ 23% &amp; Delivery.</b></p> <p style="color: red;"><b>Private clients may be entitled to reclaim VAT - please enquire</b></p>					SUB-TOTAL		
					<b>TOTAL PURCHASE PRICE</b>		

## 1. Clinical Appropriateness

### Indications for Use:

1. Neurological conditions including but not limited to stroke, brain injury, spinal cord injury, MS, GBS, ALS, and muscular dystrophy.
2. Brachial plexus injury
3. Radial nerve injury
4. Individuals exhibiting wrist drop and/or lacking finger extension.



### Contraindications:

1. Not for use with clients exhibiting increased tone/spasticity.
2. Not for use with severe contractures or joint deformities in the fingers.
3. Not for use over open wounds and infected areas.
4. Not for use with severe oedema.

## 2. How to Measure

To maximize comfort and fit, there are (4) measurements required (2 for the Wrist Splint and 2 for the Glove Liner). The four measurements will indicate the best possible size recommendation. Once you have determined the below measurements, use the chart on page 1 to determine the appropriate recommended size.

### STEP 1: WRIST SPLINT MEASUREMENTS:




Refer to the “Wrist Splint Sizing” chart on page 1 for final measurement.

A.) **Wrist Circumference:** measure wrist circumference just proximal to head of ulna: \_\_\_\_\_

B.) **Dorsal Hand Width:** measure dorsal hand width just proximal to MCPs (do not wrap tape down lateral or medial side of hand): \_\_\_\_\_

### STEP 2: GLOVE LINER MEASUREMENTS:



Refer to the “Glove Liner Sizing” chart on page 1 for final measurement.

C.) **Hand Circumference:** measure hand circumference by wrapping tape circumferentially around hand: \_\_\_\_\_

D.) **3rd Floor Length:** measure from tip of 3<sup>rd</sup> digit to the most distal wrist crease: \_\_\_\_\_