




Memo-Pad Cushion

Comfort & Care

The **Memo-Pad Cushion** has been designed for users considered to be at **'High Risk'** of ulcer development. The temperature sensitive foam conforms to the user's shape, creating a greater body contact area, thus effective pressure redistribution is achieved. The cover is both multi-stretch and vapour permeable. The high frequency welded seams and a fully concealed zip reduces the possibility of fluid ingress and aids in infection control.



Key Features

- High specification foam contributes to an effective pressure ulcer prevention plan.
- Smooth multi-stretch cover allows the patient to glide over the surface of the cushion.
- The waterproof, vapour permeable cover can be easily wiped clean.
- Cover is impregnated with an anti-fresh biocide to help prevent the build up of bacteria.
- High frequency welded seams reduce the risk of fluid ingress.
- Full length welded flap cover can be removed and laundered at 95 degrees Celsius.

Sizes

APH012 : 43cm x 43cm x 7.6 cm
APH013 : 51cm x 51cm x 7.6 cm



Specifications

Product Code: APH012 / APH013
Max. Weight Limit: 24 stone / 152.7 Kg
Risk Factor: High Risk
Warranty: 2 years
Cushion Weight: 3 Kg
Safety Standards: BS EN 597-1 Cig, BS EN 597-2 Match, BS7177: 1996 Crib Source 5

